

**MATCON**<sup>®</sup>

POWDERS, HANDLED.

# BLENDING SOLUTIONS FOR PHARMACEUTICALS

BLENDING ANY BATCH SIZE.  
EFFICIENTLY.

---


[www.matconibc.com](http://www.matconibc.com)

# The Pharmaceutical Blender for fast, efficient blending

**Improve blending efficiency by using just the one blender.**

Because the Matcon blender can accommodate different sizes of IBC, there is no need to invest in several different mixers to match batch sizes.

As the container is the blending vessel there is no need to clean between batches either.



Comes with interlocked safety cage

## THE MATCON BLENDING ADVANTAGE

System is always blending – filling, emptying and cleaning are done offline


In-bin blending for instant formulation changes without risk of contamination

Wide range of batch sizes give flexibility

Range of IBC sizes to suit batch requirements



The IBC Blending module enables the mixing of any recipe, at any time.



**Eliminate the time-consuming cleaning associated with traditional static mixers.**

**Rapid batch changeovers** as there is no product contact with the blender.

**No cross-contamination risk** as the IBC is the blending vessel and is fully contained.

**High OEE rates** as IBCs are filled, emptied and cleaned offline, reducing downtime.

**Flexible batch sizes** can be blended, as one blender accommodates a multiple of different IBC sizes.

STRUCTURE SIZES ARE DEPENDENT UPON IBC SIZE

# MATCON®

POWDERS, HANDLED.

We specialise in providing complete materials handling solutions in the Food, Nutrition, Pharmaceutical and Chemical sectors.

We're not just an IBC provider, but a true partner helping you accomplish the right system for your needs.



THE POWDER HANDLING EXPERTS

[www.matconibc.com](http://www.matconibc.com)

MATCON, Bramley Drive, Vale Park West, Evesham, Worcestershire, WR11 1JH, UK

Tel: +44 (0) 1386 769 000 Fax: +44 (0) 870 191 1116

